SARTURIUS

Sartoflow®4500

Single-use Tangential Flow Filtration System



Product Information

The Sartoflow® 4500 single-use tangential flow filtration (SUTFF) system offers all of the controls and instrumentation for running fully automated UF | DF processes for batch volumes of 200 to 2,000 L (1 to 15 kg). The system is compact, has low hold-up volume, and the capacity to run cassette surface areas from 3 to 10 m² and hollow fiber modules of 2.69 and 5 m² filtration area.

Benefits

- Designed for efficiency in SUTFF applications
- High level of operation-based automation
- User-friendly and quick set-up time

Relevant Processes

- Vaccines
- Monoclonal antibodies
- Recombinant proteins

Relevant Process Steps

- Ultrafiltration
- Diafiltration

Related Process Scale


Large scale manufacturing

System Concept

The SUTFF system offers a combination of hardware, control system and single-use flow kits designed to work together perfectly. A smart system design ensures ease of use when running fully automated TFF process sequences with complete single-use flow kits incorporating all critical sensor and instrumentation technologies.

- The SUTFF system facilitates the automated processing of batch volumes from 200 to 2,000 L of monoclonal antibodies, recombinant proteins and vaccines processes.
- The single-use sensors provide the means to fully monitor and control the process.

Control System



Main screen with indication of all process signals as well as system | component status

The automation concept is based on the Process Step Editor (PSE), which allows the user to configure automated process sequences (operations) in a very simple and user-friendly way. This is an intuitive control system that is ideally suited for use in the SUTFF system and has an operation editor that allows easy generation and configuration of end-user-specific processes.

- The user can select the transition points that move the process from one step to the next by selecting one or more process parameters (with an "and | or" logic) such as conductivity, retentate volume, permeate volume, bag weight.
- Further flexibility is added as the system allows the configuration and use of loop and jump functions, making selective step repetitions feasible. Once selected by the operator, batch reports are automatically generated for each step at the end of a batch

Quality Standards

System

A stringent approach to quality of purchased and manufactured components is maintained. The system hardware is designed and built to well recognized industry standards, including:


- Good automated manufacturing practice (GAMP) current version
- Code of Federal Regulations 21 Good manufacturing practice (GMP) for finished pharmaceuticals
- 2004/108/EC Software complies with CFR21Part11 for electronic records, signature and audit trail
- The quality management system is ISO 9001 compliant.

Single-use Flow Kits for Flat Sheet Cassettes

- Manufactured in a controlled environment (Class 7 in operation according to ISO 14644)
- Feature components tested for biocompatibility and certified according to USP<88>
- Components are compliant with EMA410/01 Rev3 requirements with regards to TSE | BSE
- Manufactured according to a validated assembly process
- Supplied double-bagged and irradiated with a minimum dose of 25 kGy
- The quality management system is ISO 9001 compliant



Sartoflow® 4500 Cassette Holder 3 to 10 m²



Single-use Flow Kits

Technical Specifications

System Specifications				
Filtration area	Flat sheet cassettes 3 – 10 m²			
Hollow fiber modules	2.69 and 5 m ²			
Recirculation pump	Quattroflow™* 5050SU 4-piston diaphragm pump 100 – 4,500 L/hr			
Fed batch buffer pump	Quattroflow™ 1200SU 4-piston diaphragm pump 20 to 900 L/hr			
Over pressure protection	Yes			
Maximum pressure	4 bar 58 psi			
Volume in 1" feed and retentate lines without casssettes	2.8 L			
Non recoverable volume	<100 mL in QUT (system flow kit , not including the retentate vessel)			
Operating temperature range	4 - 40 °C			
Weight of empty system	420 kg			
Pneumatic supply	6 bar 90 psi			
System dimensions (length × depth × height)	EMEA Asia: 1,920 × 1,000 × 1,718 mm (nominal)	US: 1,980 × 1,120 × 1,730 mm (nominal)		
Power	EMEA Asia: 400 VAC, 3-phase 50 Hz	US: 480 VAC, 3-phase 60 Hz, 11.1 FLA		
Component Specification				
Recirculation pump	100 - 4,500 L/hr			
Fed batch buffer pump	20 - 900 L/hr			
Pneumatic pinch valves	For use with 1" ID tubing (feed and retentate	For use with 1″ ID tubing (feed and retentate)		
Control valve	Single-use diaphragm valve body with posit	tioner		
Recovery valve	Single-use diaphragm valve body			
Retentate flow measurement	Mathematically derived from pump (factory o	calibration range: 900 – 1,800 L/hr accuracy ± 10% of reading)		
Permeate flow measurement	± 6% of full scale (1,000 L/hr)			
Pressure measurement	± 0.15 bar ± 2 psi of reading			
	± 3% from 0 - 0.41 bar ± 5% from 0.41 - 2.07 bar ± 7% from 2.07 - 4.14 bar			
Temperature measurement	2 - 50°C ± 0.5 °C			
Conductivity	Retentate: 0.0 to 150.0 mS/cm Permeate: 0.0 to 100.0 µS/cm and 0.1 to 15 ± 0.5 mS/cm from 0 - 10 mS/cm ± 2 mS/cm from 10 - 100 mS/cm ± 4 mS/cm from 100 - 150 mS/cm (dependent of process and ambient tempe			

^{*}Quattroflow $^{\text{TM}}$ is a trademark of Dover Company.

Component Specification			
UV	0 - 3 AU ± 0.06 AU		
рН	0 - 12 pH ± 0.2 pH		
Level sensor	Presence absence liquid on the produc	Presence absence liquid on the product inlet line	
Airflow	EMEA Asia:	US:	
	0 – 12,500 sccm	0-20,000 sccm	
	± 0.18% of FS (< 20% FS)	± 0.18% of FS (< 20% FS)	
	± 0.9% of S.P. (20 – 100% FS)	± 0.9% of S.P. (20 – 100% FS)	
Mass measurement	Appropriate size load cell will be defined	Appropriate size load cell will be defined for each process	

Components in Fluid Contact	Material of Construction
Tubing	Platinum cured silicone, reinforced in high pressure sections
Product recirculation pump	Machined polypropylene head, EPDM valves and O-rings, Santoprene™* diaphragm
Diafiltration fed batch pump	Molded polypropylene head, polypropylene ports, EPDM valves and O-rings, Santoprene™ diaphragm
Connectors and fittings	Polypropylene
Single sensors: Pressure, UV Flow Conductivity (incl. temperature)	PolysulfonePVDF with ruby bearingPolysulfone with stainless steel electrode
Multi sensor cell: Temperature conductivity UV	PPSU (electrodes: 316 L, gaskets: EPDM, UV window: quartz)
Sanitary fittings	Polypropylene, gaskets silicone
TC gaskets	Silicone
Diaphragm valves	Body polypropylene, diaphragm TPE
Distribution plate	Polypropylene

^{*}Santoprene $^{\text{\tiny{TM}}}$ is a trademark of Exxon Mobil Corporation.

Ordering Information

Preconfigured Sartoflow® 4500 System

	EMEA and APAC	US
Confoguration code	LGRTFF4500SPPUTX	WH3SUTFFSTD-1IN
System components	Quattroflow™ 5050SU recirculation p	ump (600 – 4500 L/hr) and Quattroflow™ 1200SU buffer pump
	Pneumatic pinch valves: SU diaphragm control valve SU diaphragm recovery valve	
	Single sensors: Pressure, UV, temperature, flow, cond	luctivity
	Multi sensor cell: Temperature conductivity pH U Flow, TMP and bag weight controlle	

Single-Use Flow Kits for Flat Sheet TFF Cassettes

Status 2022 Q1

Description	Sensors	Order Codes EMEA and Asia	Order Codes US
Sartoflow® STD 4500 1" - Feed	Pressure	7444-0653J	8444-0687K
Sartoflow® STD 4500 1" - Diafiltration	Pressure	7444-0653G	8444-0734H
Sartoflow® STD 4500 1" - Buffer & WFI	-	7444-0653H	8444-0687J
Sartoflow® STD 4500 1" - Retentate	Pressure, conductivity, temp.	7444-0653K	8444-0687L
Sartoflow® STD 4500 ½" - Permeate	Pressure, flow, conductivity, UV, pH, temp.	7444-0653L	8444-0734J
Sartoflow® STD 4500 - Pump Head	-	7444-0653M	8444-0734K

Germany

August-Spindler-Strasse 11 37079 Goettingen Phone +49 551 308 0

R For further contacts, visit www.sartorius.com

USA

Sartorius Stedim Biotech GmbH Sartorius Stedim North America Inc. 565 Johnson Avenue Bohemia, NY 11716 Toll-Free +1 800 368 7178